

McStas and McStasScript

3rd November, 2019

Mads Bertelsen, Thomas Holm Rasmussen
European Spallation Source




**EUROPEAN
SPALLATION
SOURCE**






McStas - Overview

- Name: Monte Carlo Simulation of Triple Axis Spectrometer
- Now simulates all kinds of neutron scattering instrumentation
- From neutron source to detector
- Core developers at DTU / ESS
- Open source community project, many contributions




McStas - Overview

- 3D model and physics
- Construct instrument from pool of components
- Each component describe spatially separate part of the instrument




McStas - Components

- **Sources**
 - From file
 - Analytical
- **Optics**
 - Guides
 - Choppers
 - Monochromator
- **Sample**
 - Crystalline
 - Inelastic
 - Large scale structures
- **Monitors**




McStas – Technical

- **Instrument file (meta language)**
 1. Instrument definition
 2. Declare section
 3. Initialize section
 4. Trace section
 5. Finally section
- **Instrument to data with one utility**
 - mcrun on command line
 - mcgui graphical user interface




McStasScript

- **Generate instrument file from python API**
- **Advantages:**
 - for loops
 - Functions
 - Flexible code placement



McStasScript

- **Generate instrument file from python API**
- **Advantages:**
 - for loops
 - Functions
 - Flexible code placement
- **Read component library**
- **Execute simulation**
- **Load data**



McStasScript

- **Installation**

Python –m pip install McStasScript

- **Configuration**

- Set path to mcrun utility

- Set path to components

```
In [2]: from mcstasscript.interface import instr, plotter, functions
```

```
In [3]: my_config = functions.Configurator()  
my_config.set_line_length(72)  
my_config.set_mcrun_path("/Applications/McStas-2.5.app/Contents/Resources/mcstas/2.5/bin/")  
my_config.set_mcstas_path("/Applications/McStas-2.5.app/Contents/Resources/mcstas/2.5/")
```



McStasScript

- Create a new instrument object

```
In [3]: demo = instr.McStas_instr("presentation_demo")
```



McStasScript

- Add instrument parameter
- Add a component
- Print component state

```
In [4]: demo.add_parameter("wavelength", value=4.0, comment="[AA] Wavelength simulated")
demo.show_parameters()
```

```
wavelength = 4.0 // [AA] Wavelength simulated
```

```
In [5]: src = demo.add_component("source", "Source_simple")
src.xwidth = 0.11
src.yheight = 0.11
src.focus_xw = 0.05
src.focus_yh = 0.05
src.dist = 2.0
src.lambda0 = "wavelength"
src.dlambda = "0.1*wavelength"
src.flux = 1E13
```

```
In [6]: src.print_long()
```

```
COMPONENT source = Source_simple
  yheight = 0.11 [m]
  xwidth = 0.11 [m]
  dist = 2.0 [m]
  focus_xw = 0.05 [m]
  focus_yh = 0.05 [m]
  lambda0 = wavelength [AA]
  dlambda = 0.1*wavelength [AA]
  flux = 10000000000000.0 [1/(s*cm**2*st*energy unit)]
AT [0, 0, 0] ABSOLUTE
```



McStasScript

- Include code in the simulation, declare and initialize

```
In [7]: monochromator_Q = 1.8734
demo.add_declare_var("double", "mono_Q", value=monochromator_Q)
demo.add_declare_var("double", "wavevector")
demo.append_initialize("wavevector = 2*PI/wavelength;")

demo.add_declare_var("double", "mono_rotation")
demo.append_initialize("mono_rotation = asin(mono_Q/(2.0*wavevector))*RAD2DEG;")
```



McStasScript

- Include code in the simulation, declare and initialize
- Use in component

```
In [7]: monochromator_Q = 1.8734
demo.add_declare_var("double", "mono_Q", value=monochromator_Q)
demo.add_declare_var("double", "wavevector")
demo.append_initialize("wavevector = 2*PI/wavelength;")

demo.add_declare_var("double", "mono_rotation")
demo.append_initialize("mono_rotation = asin(mono_Q/(2.0*wavevector))*RAD2DEG;")
```

```
In [8]: mono = demo.add_component("mono", "Monochromator_flat")
```

```
In [9]: mono.zwidth = 0.05
mono.yheight = 0.05
mono.Q = monochromator_Q
mono.set_AT([0, 0, 2], RELATIVE="source")
mono.set_ROTATED([0, "mono_rotation", 0], RELATIVE="source")
```



McStasScript

- Add a monitor to collect data

```
In [10]: beam_direction = demo.add_component("beam_dir", "Arm", AT_RELATIVE="mono")
         beam_direction.set_ROTATED([0, "mono_rotation", 0], RELATIVE="mono")

         detector = demo.add_component("detector", "PSD_monitor")
         detector.xwidth = 0.1
         detector.yheight = 0.1
         detector.filename = "\"data.dat\""
         detector.nx = 100
         detector.ny = 100
         detector.set_AT([0,0,2], RELATIVE="beam_dir")
```



McStasScript

- Run the simulation from python
- Data returned as simple object with numpy data arrays

```
In [16]: data = demo.run_full_instrument(ncount=1E7, foldername="demo", increment_folder_name=True,  
                                         parameters={"wavelength" : 3.5})
```

```
INFO: Using directory: "demo_0"
```

```
INFO: Regenerating c-file: presentation_demo.c
```

```
CFLAGS=
```

```
INFO: Recompiling: ./presentation_demo.out
```

```
INFO: ===
```

```
INFO: Placing instr file copy presentation_demo.instr in dataset demo_0
```


```
Detector: detector_I=8.35572e+09 detector_ERR=2.10199e+07 detector_N=158020 "data.dat"
```



McStasScript

```
In [12]: plotter.make_plot(data)
```

```
number of elements in data list = 1  
Plotting data with name detector
```




McStasScript

```
In [13]: functions.name_plot_options("detector", data, log=True)
```

```
In [14]: plotter.make_plot(data)
```

number of elements in data list = 1
Plotting data with name detector



Plot preferences included in data object



McStasScript

```
In [1]: from mcstasscript.interface import instr, plotter, functions, reader
```

```
In [2]: InstrReader = reader.McStas_file("presentation_demo.instr")
```

```
In [3]: InstrReader.write_python_file("generated_presentation.py", force=True)
```

```
In [4]: demo = instr.McStas_instr("generated_presentation_demo")
InstrReader.add_to_instr(demo)
```

```
In [5]: demo.print_components()
```

source	Source_simple	AT	(0, 0, 0)	ABSOLUTE	
		ROTATED	(0, 0, 0)	ABSOLUTE	
mono	Monochromator_flat	AT	(0, 0, 2)	RELATIVE	source
		ROTATED	(0, mono_rotation, 0)	RELATIVE	source
beam_dir	Arm	AT	(0, 0, 0)	RELATIVE	mono
		ROTATED	(0, mono_rotation, 0)	RELATIVE	mono
detector	PSD_monitor	AT	(0, 0, 2)	RELATIVE	beam_dir
		ROTATED	(0, 0, 0)	RELATIVE	beam_dir



McStasScript

- Python API as an alternative to standard McStas workflow
- Still work to do on API, may rename some methods / functions / keywords
- Improvements on plotting and data handling
- Review pending from co-workers
- Continuous integration will be set up



PaNOSC WP5 deliverables

- **openPMD**
 - Has been postponed slightly in favor of building an API
 - Have example dataset and extension document
 - Need to make McStas in/out components
- **Documented API**
 - Have an API, will need further improvements
 - Both docstrings and pdf documentation
 - 199 unit tests of core functionality, lacks test of plotting / instrument reader
 - Wish to add continuous integration
- **Release of examples in Jupyter Notebook**
 - Can generate python scripts for McStas instruments (Hundreds available)
 - Can create more examples that emphasises python advantages



Thank you

Mads.Bertelsen@esss.se

